



Medical Stool

Superior infection control features

Designed from top to bottom for the medical profession.

The **healthCentric Medical Stool Series** offers superior infection control features. The ring height adjustment activator ensures ease of adjustability from every angle. Additional base and lift options provide the ultimate flexibility for operating rooms, examination areas, nursing stations or wherever mobile seating is required.

Superior Infection Control

Proprietary to healthCentric, IC+™ Upholstery Solution* is a breakthrough technology specially formulated for the healthcare industry. Providing durability where it matters most, IC+ has a superior abrasion rating of 1,000,000 double rubs while guaranteed not to crack, puncture or tear for 10 years.

- + Puncture and tear resistant
- + Sealed and seamless
- + Hospital-grade cleaning agent tough†
- + Non-porous
- + Guaranteed not to crack
- + Bed bug proof
- + PVC- and latex-free
- + GREENGUARD® certified
- + 10 year warranty
- + CAL 133 compliant



Upholstered / Vinyl Seat Style



IC+ Seat Style



Features

Seat

- 16" Diameter upholstered seat (fabric, vinyl, or IC+)
- 13" Polyurethane seat (PU)
- Height adjustment
- Easy-to-use treadle lever height adjustment activator

Pneumatic lift

- 200 mm Knurled dual height adjustment system
- Easy-to-use ring type height adjustment activator
- Adjustment between 22.75" - 29.75"
- Self-skinning Polyurethane (PU) seat reduces seat height by 2".

Base

- 22" Chrome 22 Base
- 18" Oval tube footring built in for durability
- Dual wheel urethane casters

Options

- IC+ Under-Seat Seal
- 2 1/4" Glides
- Lock-On Urethane Casters
- Lock-Off Urethane Casters
- Single Wheel Urethane Casters
- CAL 133

Warranty

LIFETIME°

- Original purchaser
- Defects and functional failure of materials up to a weight of 280 lbs
- With the following exceptions:

10 Years°

- Foam, fabric
- IC+ against cracking, punctures and tears

5 Years°

- Foot treadle lever mechanism

° Includes labor costs.

* IC+ Upholstery Solution patent pending.

† When used in the quantities and manners specified by the cleaning agents' labeled cleaning instructions.

